


# HOW TO USE YOUR NEW FEEDING PUMP

Your new pump is there to make things a bit easier for you. Remember, using your pump shouldn't mean you have to miss out on anything you like to do. Even if you are not allowed to eat normal food, you can still have **lots** of fun!



## WHAT WORDS MEAN

Before you start using your pump, learn about some strange new words

|                          |  |
|--------------------------|--|
| Volume, Flow rate, Speed | Tells you how <b>fast</b> your pump is running   |
| Volume, Dose             | Tells you how <b>much</b> feed is given  |
| Giving set               | This is the tubing that you connect to your pump and tube feed. It goes through the pump and connects to your feeding tube |
| Continuous feeding       | A feed given over a long time, like lots of meals joined into one  |
| Bolus feeding            | A feed given over a short time, like one meal  |



YOUR PUMP IS  
SMALL, YOU CAN  
TAKE IT WITH YOU  
WHEREVER YOU GO



IF YOU DON'T HAVE A  
BACKPACK CONTACT  
YOUR NUTRICIA  
HOMEWARD NURSE

# HOW TO SET UP YOUR PUMP

Your new pump is quite easy to use — these steps will teach you how to use your pump. Please ask an adult if you need help at any time.

## STEP 1


Open the feed bottle; take the giving set and screw it tightly onto the bottle, so the spike goes through the silver foil; hang the feed up



## STEP 2

Open the pump door by pressing on the door clip (A) and lift the door (B) at the same time

Put the loop part of the giving set around the big wheel, pulling slightly on the giving set to clip the plastic part into place; close the pump door




## STEP 3



Press the **ON/OFF** button and hold it down, count slowly to three; the pump is on when you hear a beep



Press the **CLR** button to clear the volume memory from your last feed



Press the **FILL SET** button and hold it down, count slowly to three; the pump will start and the whole tube will fill up with the feed. You can stop the tube filling up with feed by pushing the FILL SET button again

**YOU ARE NOW READY TO SET UP YOUR PERSONAL FEED PROGRAMME**

## HOW TO SET UP A CONTINUOUS FEED

- Check that your pump is switched on
- Press the **ml/h** button, then use the **+/-** buttons to set the 'rate' (or speed) your pump will run at — this is the rate your dietitian or nurse has asked you to use
- Press the **START/STOP** button to start feeding

**IF YOU NEED TO STOP YOUR PUMP**, press the **START/STOP** button once; the pump will stop when you hear three beeps

**TO SWITCH YOUR PUMP OFF** when the feed has finished, press the **ON/OFF** button and hold it down, counting slowly to three until you hear a beep — the pump is now off

## HOW TO SET UP A DOSE FEED

- After setting the rate, press the **DOSE=VOL** button
- Use the **+/-** buttons to set the dose (or total volume to be delivered)
- Press the **START/STOP** button to start feeding

Read through the steps again when you need to and don't forget that your dietitian or nurse will be happy to answer any questions you might have.

Happy feeding, and don't let your pump stop you from having fun!



Nutricia Advanced Medical Nutrition  
White Horse Business Park, Trowbridge, Wiltshire BA14 0XQ

SCC3781-12/19