

NASO-JEJUNAL (NJ) TUBE

THE FOLLOWING INFORMATION IS AIMED AT SUPPORTING YOU AND YOUR CHILD ON HOW TO CARE FOR THEIR NJ TUBE.

IMPORTANT NOTE: It is important that you contact your Nutricia Homeward Nurse if you need additional training on the Flocare pump (if applicable) or tube feeding.



DAILY CARE OF YOUR CHILD'S NJ TUBE

It is important to keep the tube in good condition to avoid unnecessary replacement of the tube. Regular care and flushing of the tube will help to prevent the tube becoming blocked.

- Wash hands before and after handling the tube
- Confirm the position of the NJ tube by checking the cm graduation marking at the nostril remains the same as noted after initial placement. Position must be checked before administration of any water, enteral nutrition or medication. Never administer anything down the tube before confirming the position
- Flush the tube as per training provided by Nutricia Homeward Nurse using a 60ml enteral syringe with water (type and volume as recommended for jejunal use by the managing healthcare professional) to prevent the tube from blocking. Water should be administered before and after enteral nutrition or medication administration or every 4-6 hours if feeding is not in progress (except during the night) or as recommended by the managing healthcare professional
- Care of the nose: fixation tape should be changed weekly, but more frequently if it loosens or becomes dirty. When changing the tape, clean skin thoroughly and try to alter the position of the tape from where it was last fixed. If the skin or nostril becomes sore or red contact a healthcare professional for instructions on how to proceed. Avoid the use of creams or powders as they can damage the tube
- Never attempt to replace the NJ tube as this type of tube can only be replaced by a healthcare professional in the hospital
- Check the naso-jejunal tube daily to ensure it is intact and that the tube is not kinked or blocked and remains securely fixed to the face.

HOW DO I CHECK THE POSITION OF MY CHILD'S NJ TUBE?

It is important to check the tube position prior to administering anything via the feeding tube, after a coughing fit or following a vomiting episode.

- Confirm the position of the tube by checking the cm graduation marking at the nostril remains the same as noted after initial placement in the hospital. The position of the tube must be checked before administration of any water, enteral nutrition or medication
- If there is any doubt about the position of the tube, do not use the tube and contact the managing healthcare professional immediately for instructions.

IMPORTANT NOTES:

You are advised to contact your managing healthcare professional immediately if there is any reason that medication or enteral nutrition needs to be delayed or omitted due to no available enteral or oral route.

It is important your child remains upright at a 45° angle during enteral feeding and for at least one hour after enteral nutrition has been administered.

TROUBLESHOOTING NASO-JEJUNAL TUBES

MY CHILD'S NJ TUBE HAS BECOME BLOCKED

If there is resistance when flushing the naso-jejunal tube, do not force water into the tube:

- Confirm position of the NJ tube before administering anything down the tube
- Using a gentle pull and push technique, flush the tube using a 60ml enteral syringe with warm water suitable for jejunal use (type and volume as recommended by the managing healthcare professional). **Do not use acidic solutions such as fruit juices or cola as they can curdle the enteral nutrition product**
- If unable to clear the blockage contact the managing healthcare professional or Nutricia Homeward Nurse for instructions on what to do.

MY CHILD'S NJ TUBE HAS COME OUT

It is important that your child has a new NJ tube placed in time for their next enteral nutrition or medication administration.

A NJ tube must always be replaced in a hospital setting therefore do not attempt to replace a new tube at home.

Follow the steps below:

- Remain calm
- Contact the managing healthcare professional and explain that the naso-jejunal tube has come out
- If the managing healthcare professional is unavailable contact your GP or local hospital emergency department to let them know the NJ tube will need to be replaced. This will give the department time to find the medical notes and to ensure a member of staff is available to replace the NJ tube
- Take the tube that has fallen out with you in a clean plastic bag so that the staff can identify which type of tube is required.

The Nutricia Homeward Nursing Service advice literature is provided as guidance following Nutricia Homeward Nursing Service training, and should not be used as a substitute for medical advice. Always contact your healthcare professional prior to making any changes to your tube feeding regimen.

Please contact the Nutricia Homeward Nurse or managing healthcare professional for the most up to date version every year. Provided on behalf of Nutricia Ltd at the request of the NHS managing healthcare professional.

IMPORTANT: Always seek urgent medical advice if there are any signs of abdominal pain, chest pain or breathing difficulties.